

Arizona Data Qualifier Rev. 4.0 ADEQ Accpetance Matrix

A	B	C	D	E	H	K	L	M	N	Q	R	S	T	V
A1	B1	C1	D1	E1	H1	K1	L1	M1	N1	Q1	R1	S1	T1*	V1
A2	B2		D2	E2	H2	K2	L2	M2	N2	Q2*	R2		T2*	V2
A3	B3	C3	D3	E3	H3		L3	M3		Q3*		S3	T3*	V3
A4	B4	C4	D4	E4	H4		L4	M4	N4	Q4	R4	S4	T4	
A5	B5	C5	D5	E5	H5	K5	L5	M5	N5	Q5	R5	S5	T5	V5
A6	B6		D6	E6	H6	K6		M6	N6	Q6	R6	S6	T6	V6
A7*	B7		D7*	E7		K7		M7	N7	Q7	R7	S7	T7	
A8		C8		E8		K8				Q8	R8	S8		
						K9				Q9	R9			V9
						K10				Q10		S10		
						K11				Q11	R11	S11		
											R12	S12		
											R13			

Acceptable/Lab Discretion	<p>A7* - Must transfer to sterile container to homogenize D7* - MRL must < MCL</p> <p>Q2* - not acceptable for VOCs or TTHMs Q3* - only acceptable if preserved within 24 hrs of sampling</p> <p>T1, T2, T3* - ADHS will consider approving new methods upon request from licensed lab.</p>
Acceptable if ND	
Acceptable if > MCL	
Not Acceptable	
	Qualifier Deleted